

Y	ELECTRONICS			
	AES System Upgrades			
Software Upgrades	Description		Notes	
	SmartSoft Auger for Windows 10 PC upgrade for 680,690, 700, 700Xi	NA	Requires Hybrid Electronics Upgrade	Sales
	SmartSoft Auger for Windows 10 PC upgrade 710 and Hybrid	Price book		Sales
Hardware erformance	680,690, 700, 700Xi, 710 Transfer Vessel	Price book		Sales
	Hybrid Electronics with optional Hero mode and intro camera	Price book	680 and above	Sales
	710 BSE Upgrade	Special		Sales
	710 FIB Upgrade	Price book	Requires large enclosure	Sales
Harc Perfol	710 EDS Upgrade	Price book	Requires large enclosure	Sales
Pe H	710 EBSD Upgrade	Price book	Maybe possible for 700Xi depending on the chamber configuration	Sales
	680,690 EDS Upgrade	Price book		Sales
	XPS System Upgrades			Requested to:
	Description	/	Notes	
	VersaProbe Hybrid, II and III SmartSoft for Windows 10 PC Upgrade	Price book	Notes	Sales
es	Quantera Compass to SmartSoft Upgrade	NA	Requires Hybrid Electronics upgrade	Sales
wa	Quantera Compass to SmartSoft Opgrade Quantes, Quantera Hybrid and II SmartSoft for Windows 10 PC Upgrade	Price book	Requires Trybrid Electronics apgrade	Sales
Software Upgrades	Image registration for VersaProbe Hybrid, II, III, Quantera Hybrid and II	Price book	Requires SmartSoft for Windows 10 or 11	Sales
S D	StrataPHI Thin Film Structure Analysis Software	Price book	Requires SmartSort for Windows 10 of 11	Sales
	Statut III IIII Statutare / Mary 55 Software	THEE BOOK		34163
	VersaProbe I			
	VersaProbe to VersaProbe Hybrid	Price book	Optional fixed Leak	Sales
	Transfer Vessel	Special	May require an adapter (Special)	Sales
	VersaProbe Hybrid, II and III:	Special	iviay require arradapter (Speciar)	Sales
	VP III- AD optics	Price book		Sales
e e	UV Source	Price book	New UV source does not require additional pumps	Sales
var	Dual Anode	Price book	New ov source does not require additional pamps	Sales
Hardware Performance	LEIPS	Price book		Sales
Hal	REELS	Price book		Sales
مّ	Scanning Auger	Price book	Requires Mu metal chamber (not possible for VP Hybrid with 06-350)	Sales
	Hot/Cold Stage	Price book	,	Sales
	Heating only Stage	Price book		Sales
	4 Contact H/C Stage	Special		Sales
	H/C intro	Special	Requires 60mm intro	Sales
	60mm intro	Special	Strongly recommended for 4 Contact H/C stage	Sales
	Intro parking	Special	Requires 60mm intro	Sales
	Intro dual turbo	Special	Recommended for intro parking	Sales
	Main Chamber Turbo	Price book	Required with GCIB for system with ion pump in analysis chamber	Sales
	20 keV GCIB	Price book		Sales
	20 keV C60	Price book		Sales
a	Customized preparation chamber	Special		Sales
vare nance	GCIB cluster measurement tool	Price book	Only available for GCIB with 2 turbos	Sales
_ ×a ma	Transfer Vessel	Price book	25mm and 60mm	Sales

Page 1 of 2 January, 2024



Hard Perfor	Quantera			
	Quantera to Quantera Hybrid Upgrade	Price book		Sales
	Quantera Hybrid and II			
	Auxiliary chamber	Special		Sales
	Secondary intro with Transfer Vessel (25 and 60mm) and optional camera	Special		Sales
	H/C Intro and Stage	Special	Also with contacts for in-situ electrochemical experiments	Sales
	Dual Turbo	Special		Sales
	Glove box interface for standard intro	Special		Sales

Special

4 contacts ambient stage

For in-situ electrochemical experiments

	SIMS System Upgrades							
	Description		Notes					
	ADEPT 1010 Windows 10 PC Update	Price book		Service				
es e	ADEPT 1010 Hybrid Upgrade	Price book		Sales				
Software Upgrades	ADEPT microbeam Upgrade	Special		Sales				
	nanoTOF I Win 10 PC Update	Price book		Service				
	nanoTOF II SmartSoft-TOF for Windows 10 PC Upgrade	Price book		Sales				
	Windows 10 TOF-DR Data reduction software (nanoTOF I & II)	Price book		Sales				
	nanoTOF I							
	nanoTOF I to nanoTOF II	Special		Sales				
	nanoTOF II							
9	MS/MS (Ready)	Special		Sales				
ano	MS/MS (Non-Ready)	NA		Sales				
Hardware Performance	Pulsed option for monoatomic ion gun	Price book		Sales				
rfo	H/C stage and intro	Special		Sales				
Pe	Heating only stage and intro	Special		Sales				
are	Cs gun	Special		Sales				
<u>×</u>	20 keV GCIB gun	Special		Sales				
ard	20 keV C60 gun	Special		Sales				
至	FIB gun (optional spectroscopy capability)	Price book		Sales				
	Fast Rotation Module	Special		Sales				
	nanoTOF I, II Transfer Vessel Adapter	Price book		Sales				
	nanoTOF I, II Transfer Vessel	Price book		Sales				

Page 2 of 2 January, 2024

Sales